

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

9-Ethyl-3-carbazolecarboxaldehyde - 98%

Product Number:**151483**

CAS Number:

7570-45-8

MDL:

MFCD00004963

Formula:


C₁₅H₁₃NO

Formula Weight:

223.27 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

Conforms to Requirements

White to Yellow to Brown or Green

Appearance (Form)

Conforms to Requirements

Powder, Crystals, Crystalline Powder or Chunks

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 97.5 %

Specification: PRD.0.ZQ5.10000042044

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.