

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Manganese - chips, thickness <2.0 mm, 99%

Product Number: 266167
CAS Number: 7439-96-5
MDL: MFCD00011111
Formula: Mn
Formula Weight: 54.94 g/mol

Mn

TEST	Specification
Appearance (Color)	Brown to Black
Appearance (Form) Chips	Conforms
Thickness < 2.0 mm	Conforms
ICP Major Analysis Confirms Manganese Component	Confirmed
Complexometric EDTA	98.5 - 101.5 %

Specification: PRD.1.ZQ5.10000017970

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.