

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification


Product Name:
Lithium tri-tert-butoxyaluminum hydride solution - 1.0 M in THF

Product Number: 241814

MDL: MFCD00011532

Formula: C₁₂H₂₈AlLiO₃

Formula Weight: 254.27 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Colorless
Appearance (Form)	Liquid
Hydrogen Evolution	0.95 - 1.10 M

Specification: PRD.2.ZQ5.10000005682

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.