

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Heparin Cofactor II from human plasma – aqueous glycerol solution, ≥90% (SDS-PAGE)

Product Number:

H0417

Storage Temperature:

-20 °C

TEST	Specification
mg/ml	≥ 5.0
Protein concentration by UV Absorbance at A280	
units/mg Protein	≥ 700
One unit will inhibit one NIH unit of Thrombin.	
Purity by SDS Electrophoresis	≥ 90 %
Formulation	Conforms
Solution in 50% glycerol	
Tested For Infectious Agents	Tested

Specification: PRD.3.ZQ5.10000013381

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.