

## Product Specification

**Product Name:**


Bleomycin sulfate from *Streptomyces verticillus* - for fluorescence, mixture of bleomycin sulfate salts, lyophilized, powder or crystals, white to off-white

**Product Number:**

15361

**CAS Number:**

9041-93-4

**Storage Temperature:**

2 - 8 °C

TEST	Specification
Appearance (Form)	Powder or Crystals
Appearance (Color)	White to Off White
Solubility (Solvent)	Water
Solubility (Conc)	20.0 - 22.0 mg/ml
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
Copper (Cu)	≤ 0.1 %
Loss on Drying	≤ 6.0 %
Activity	1.5 - 2.0
U/mg	
Purity (HPLC)	≥ 90 %
(A2+ B2, anhydrous basis)	
Recommended Retest Period	-----
7 Years	

Specification: PRE.1.ZQ5.10000013080

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

