

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)

Email USA: [techserv@sial.com](mailto:techserv@sial.com)

Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

**Product Specification**

Product Name:

1-Methoxy-2-propanol ≥99.5%

Product Number:

**484407**

CAS Number:

107-98-2

MDL:


MFCD00004537

Formula:

C<sub>4</sub>H<sub>10</sub>O<sub>2</sub>

Formula Weight:

90.12 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥99.5 %
2-Methoxy-1-Propanol	≤0.49 %
Water (by Karl Fischer)	≤0.10 %
Color Test	≤10 APHA
Supplier information Product of Dow Chemical	Conforms

**Remarks:**

Specification Date : 06/21/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.