

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Nickel - foil, thickness 0.5 mm, 99.98% trace metals basis

Product Number: 357553
CAS Number: 7440-02-0
MDL: MFCD00011137
Formula: Ni
Formula Weight: 58.69 g/mol

Ni

TEST	Specification
Appearance (Color) Silver or Grey	Conforms
Appearance (Form) Foil	Conforms
X-Ray Diffraction	Conforms to Structure
Thickness 0.5 mm	Conforms
Trace Metal Analysis	≤ 250.0 ppm
Purity 99.98% Based On Trace Metals Analysis	Meets Requirements

Specification: PRD.1.ZQ5.10000008647

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.