


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ethanolamine - liquid, BioReagent, suitable for cell culture, ≥98%

Product Number: E0135
CAS Number: 141-43-5
MDL: MFCD00008183
Formula: C₂H₇NO
Formula Weight: 61.08 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Titration by HCL	≥ 98.0 %
Cell Culture Test	Pass
Recommended Retest Period	-----
4 Years	

Specification: PRD.0.ZQ5.10000018738

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.