

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Sorbitan sesquioleate - nonionic surfactant

Product Number: **S3386**  
CAS Number: 8007-43-0  
MDL: MFCD00151568



TEST	Specification
Appearance (Clarity)	Clear to Slightly Hazy
Appearance (Color)	Faint Yellow to Very Dark Yellow
Appearance (Form)	Liquid
Purity (GC)	≥ 60.0 %
Balance Primarily Palmitic, Stearic, and Linoleic Acids	
Acid Value	≤ 13.0
Hydroxyl Value	185 - 215
Saponification Value	145 - 160
Water (by Karl Fischer)	≤ 1.0 %
Brookfield Viscosity @ 25°C	1000 - 2000 cps
Recommended Retest Period 3 Years	-----

Specification: PRD.0.ZQ5.10000046761

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.