


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Calcium chloride hexahydrate - 98%

Product Number: 442909
CAS Number: 7774-34-7
MDL: MFCD00149616
Formula: CaCl₂ · 6H₂O
Formula Weight: 219.08 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form) Powder, Crystals or Granules and/or Chunks	Conforms to Requirements
Loss on Drying	49 - 50 %
Complexometric EDTA	97.5 - 102.5 %

Specification: PRD.0.ZQ5.10000023842

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.