

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-RANTES antibody produced in rat - clone 53405.111, purified immunoglobulin, lyophilized powder

Product Number:

R9400

TEST	Specification
Formulation	Conforms
Lyophilized from a solution of PBS with 5% Trehalose	
Working Dilution	Conforms
Immunocytochemistry; 8-25 ug/mL for detection in mouse splenocytes.	
Bioactivity	Conforms
measured by its ability to neutralize the biological activity of recombinant mouse RANTES. The ND50 is typically 0.1-0.5 ug/mL in the presence of 0.05 ug/mL of mouse CCL5/RANTES.	

Specification: PRD.0.ZQ5.10000047526

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.