

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:


OptiPrep™ Density Gradient Medium – used for cell and subcellular organelle isolation

Product Number:

D1556

MDL:

MFCD00867965



### TEST

### Specification

Bulk Density

1.319 - 1.321 g/ml

Endotoxin Level

≤ 1 EU/ml

Sterility

Pass

Trademark

Confirmed

OptiPrep is a trademark of Alere Technologies AS

Specification: PRD.1.ZQ5.10000033831

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.