


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
5-Methylquinoxaline - 98%

Product Number: 272310
CAS Number: 13708-12-8
MDL: MFCD00012335
Formula: C₉H₈N₂
Formula Weight: 144.17 g/mol



TEST

Specification

Appearance (Color)

Yellow to Very Dark Yellow and Light

Appearance (Form)

Brown-Yellow to Brown-Yellow and Light Brown to Brown

Infrared spectrum

Liquid or Solid

Purity (GC)

Conforms to Structure

≥ 97.5 %

Specification: PRD.0.ZQ5.10000024407

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.