

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1-Methylimidazole – ReagentPlus®, 99%

Product Number:

M50834

CAS Number:

616-47-7

MDL:


MFCD00005292

Formula:

C₄H₆N₂

Formula Weight:

82.10 g/mol



TEST

Specification

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Color Test

≤ 400 APHA

Purity (GC)

≥ 98.5 %

Specification: PRD.0.ZQ5.10000032193

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.