


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
DL-Leu-Gly-Gly

Product Number: L4127
CAS Number: 4337-37-5
MDL: MFCD00064990
Formula: C₁₀H₁₉N₃O₄
Formula Weight: 245.28 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear to Slightly Hazy
50 mg/mL, 80 % HOAc	
Infrared spectrum	Conforms to Structure
Purity (TLC)	≥ 99 %

Specification: PRD.0.ZQ5.10000103534

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

