

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Macrophage Inflammatory Protein-1 β antibody produced in goat - IgG fraction of antiserum, lyophilized powderProduct Number: **M5917**

MDL: MFCD00213001

Storage Temperature: -20 °C

TEST	Specification
Formulation	Conforms
Lyophilized from a solution in phosphate buffered saline with 5% trehalose	
Bioactivity	Conforms
ND50: 4-12 ug/mL	
Working Dilution	Conforms
by Western Blot	
Recommended working dilution is 1 ug/mL. It detects 5 ng/lane under non-reducing and 25 ng/lane under reducing conditions	

Specification: PRD.0.ZQ5.1000021155

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.