

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Neodymium(III) carbonate hydrate - 99.9% trace metals basis

Product Number:	451185	$\text{Nd}_2(\text{CO}_3)_3 \cdot x\text{H}_2\text{O}$
CAS Number:	38245-38-4	
MDL:	MFCD00150426	
Formula:	$\text{C}_3\text{Nd}_2\text{O}_9 \cdot x\text{H}_2\text{O}$	
Formula Weight:	468.51 g/mol	

TEST	Specification
Appearance (Color)	Purple
Appearance (Form)	Powder or Crystals
Complexometric EDTA % Nd	≤ 60.4 %
ICP Major Analysis Confirms Neodymium Component	Confirmed
Purity 99.9% Based On Rare Earth Analysis	Meets Requirements
Trace Rare Earth Analysis	≤ 1500.0 ppm

Specification: PRD.0.ZQ5.10000061438

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.