

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


11-Acetylmercaptoundecylphosphonic acid - 97%

Product Number:

760145

CAS Number:

304012-57-5



Formula:

C₁₃H₂₇O₄PS

Formula Weight:

310.39 g/mol

Storage Temperature:

2 - 8 °C

TEST

Specification

Appearance (Color)

White to Beige

Appearance (Form)

Powder or Crystals

Infrared spectrum

Conforms to Structure

Purity (Titration by NaOH)

96.5 - 103.5 %

GC (area %)

≥ 96.5 %

Specification: PRD.1.ZQ5.10000026809

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.