

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1,2-Bis[(2S,5S)-2,5-diethylphospholano]benzene(1,5-cyclooctadiene)rhodium(I) tetrafluoroborate

Product Number:

698385

CAS Number:


213343-64-7

Formula:

C₃₀H₄₈BF₄P₂Rh

Formula Weight:

660.36 g/mol



TEST

Specification

Appearance (Color)

Orange to Brown

Appearance (Form)

Powder

Proton NMR spectrum

Conforms to Structure

Phosphorus NMR Spectra

Conforms to Structure

Carbon

52.4 - 56.8 %

Specification: PRD.1.ZQ5.10000012616

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.