

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

4,4'-Bis[(4-bromophenyl)phenylamino]biphenyl - 97%

Product Number:

**656674**

CAS Number:

344782-48-5

MDL:


MFCD08276784

Formula:

C<sub>36</sub>H<sub>26</sub>Br<sub>2</sub>N<sub>2</sub>

Formula Weight:

646.41 g/mol

**TEST****Specification**

Appearance (Color)

White to Light Yellow and Faint Beige to Beige

Appearance (Form)

Powder or Solid or Crystals or Crystalline Chunk(s) or Chunk(s)

Infrared spectrum

Conforms to Structure

Carbon

64.6 - 69.2 %

Nitrogen

4.1 - 4.5 %

Purity (HPLC)

≥ 96.5 %

Specification: PRD.0.ZQ5.10000026815

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.