

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Erucic acid – ≥99% (capillary GC)

Product Number: **E3385**
CAS Number: 112-86-7
MDL: MFCD00063188
Formula: C₂₂H₄₂O₂
Formula Weight: 338.57 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
at 50 mg plus 1 ml of CHCl ₃	
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 99 %
Recommended Retest Period	-----
6 years	

Specification: PRD.0.ZQ5.10000088908

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.