

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-phospho-S6-Kinase (pThr³⁸⁹) antibody produced in rabbit - affinity isolated antibody, buffered aqueous glycerol solution**Product Number:** S6311

MDL: MFCD02097170

Storage Temperature: -20 °C

TEST	Specification
Formulation Supplied in 10 mM sodium HEPES, pH 7.5, containing 150 mM sodium chloride, 100 ug/ml bovine serum albumin and 50% glycerol.	Conforms
Working Dilution for immunoblotting: 1:1,000 using serum-treated NIH3T3 cell extract, or insulin-treated HeLa, COS, and C6 cell extracts.	Conforms

Specification: PRD.1.ZQ5.1000022453

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.