

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
L-15 Medium (Leibovitz) - With L-glutamine, powder, suitable for cell culture

Product Number: **L4386**

MDL: MFCD00217482

Storage Temperature: 2 - 8 °C

TEST	Specification
Origin	Non-Animal Source
Osmolality	296 - 328 mOs/kg
Galactose	5.89 - 7.20 %
Water (by Karl Fischer)	≤ 2.0 %
pH	7.8 - 8.4
Solubility (Turbidity)	Clear
Solubility (Form)	Solution
at 13.8 g/L in water	
Cell Culture Test	
Report % of control	
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Endotoxin Level	≤ 1.0 EU/ml
at 1x	

Specification: PRD.1.ZQ5.10000038321

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.