

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Ammonium lauryl sulfate solution - ~30% in H₂O

Product Number:

681806

MDL:


MFCD00050675

Formula:

C₁₂H₂₉NO₄S

Formula Weight:

283.43 g/mol



TEST

Specification

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (Titration by NaOH)

27.0 - 33.0 %

Water (by Karl Fischer)

67.0 - 73.0 %

Specification: PRD.1.ZQ5.10000054242

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.