


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Propyl gallate – powder

Product Number: **P3130**
CAS Number: 121-79-9
MDL: MFCD00002196
Formula: C₁₀H₁₂O₅
Formula Weight: 212.20 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Solubility (Color)
Solubility (Turbidity)
50 mg/mL, EtOH
Infrared spectrum
Purity (HPLC)
Recommended Retest Period
5 years

White to Off-White
Powder
Colorless to Faint Yellow
Clear

Conforms to Structure
≥ 98 %

Specification: PRD.0.ZQ5.10000030031

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.