

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Poly(3,4-ethylenedioxythiophene)-poly(styrenesulfonate) - 3.0-4.0% in H2 O, high-conductivity grade

Product Number: 655201

CAS Number: 155090-83-8 MDL: MFCD07371079

Storage Temperature: 2 - 8 ℃

TEST Specification

 Appearance (Color)
 Dark Blue

 Appearance (Form)
 Liquid

 Solid Content
 3.0 - 4.0 %

 Conductivity
 ≥ 200 S/cm

Specification: PRD.1.ZQ5.10000030456

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.