


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Methylhydroquinone - 99%

Product Number: 112968
CAS Number: 95-71-6
MDL: MFCD00002345
Formula: C₇H₈O₂
Formula Weight: 124.14 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

Off-White to Beige

Appearance (Form)

Powder to Crystalline Powder

Infrared spectrum

Conforms to Structure

Titration

98.5 - 101.5 %

% Purity With Ce(SO₄)₂

Purity (GC)

≥ 98.5 %

Specification: PRD.1.ZQ5.10000009165

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.