


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tolylene-2,4-diisocyanate - 95%

Product Number: T39853
CAS Number: 584-84-9
MDL: MFCD00002011
Formula: C₉H₆N₂O₂
Formula Weight: 174.16 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color) White (Solid Phase) or Colorless (Liquid Phase)	White or Colorless
Appearance (Form) May be Moist (Solid Phase) May be Viscuous (Liquid Phase)	Solid or Liquid
Infrared Spectrum	Conforms to Structure
Purity (GC)	≥ 94.5 %

Specification: PRD.2.ZQ5.10000035114

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.