

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Carboxypeptidase A—Agarose – ammonium sulfate suspension, ≥6 units/mL packed gel, 25 °C,
carboxypeptidase A from bovine pancreasProduct Number: **C1261**

MDL: MFCD00130718

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

units/ml	≥ 6
packed gel	
One unit will hydrolyze 1.0 micromole of Hippuryl-L-Phenylalanine per minute at pH 7.5 at 25 Deg C.	
ml suspension/ml packed gel	1.5 - 2.4
Gel content	

Specification: PRD.0.ZQ5.10000074723

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.