

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-Leukemia Inhibitory Factor antibody produced in mouse - clone 9824, purified immunoglobulin, lyophilized powder

Product Number: **L0669**

Storage Temperature: -20 °C

TEST	Specification
Appearance (Form) from a 0.2 um filtered solution in PBS with 5% trehalose	Lyophilized Powder
Endotoxin Assay	Less than 0.1 Eu per 1 ug
Elisa Capture	2 ng/ml
Neutralization of Bioactivity	0.06 to 0.2 ng/ml
Immunohistochemistry	25 ug/ml

Specification: PRD.0.ZQ5.10000087359

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.