

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Platinum – wire, diam. 0.20 mm, ≥99.99% trace metals basis

Product Number:	357359	Pt
CAS Number:	7440-06-4	
MDL:	MFCD00011179	
Formula:	Pt	
Formula Weight:	195.09 AMU	

TEST	Specification
Appearance (Color)	Conforms
Silver	
Appearance (Form)	Conforms
Wire	
Diameter	Conforms
0.20 mm	
ICP Major Analysis	Confirmed
Confirms Platinum Component	
Trace Metal Analysis	≤ 100.0 ppm
Purity	Conforms
> = 99.99% Based On Trace Metals Analysis	

Specification: PRD.0.ZQ5.10000026943

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.