

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Celery seed oil - natural, FCC, FG

Product Number: W227102
CAS Number: 8015-90-5
MDL: MFCD00146329

TEST	Specification
Appearance (Color)	Yellow to Green to Brown
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.480 - 1.490
Infrared spectrum	Conforms to Structure
Specific Gravity	0.870 - 0.910
Acid Value	≤ 4.5 ml
Saponification Value	25.00 - 65.00
Solubility (Turbidity)	Clear to Turbid
Solubility (Color)	Pass
Angular Rotation (Observed Rotation)	48 - 78 °
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period 5 Years	-----

Specification: PRD.1.ZQ5.10000038301

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.