

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-[Tris(trimethylsiloxy)silyl]propyl vinyl carbamate - ≥94%

Product Number: **795402**
CAS Number: 134072-99-4

Formula: C₁₅H₃₇NO₅Si₄
Formula Weight: 423.80 g/mol
Storage Temperature: 2 - 8 °C



| TEST | Specification |
|------|---------------|
|------|---------------|

| | |
|----------------------|---|
| Appearance (Clarity) | Clear |
| Appearance (Color) | Colorless to Yellow and Colorless to Orange |
| Appearance (Form) | Liquid or Viscous Liquid |
| Proton NMR Spectrum | Conforms to Structure |
| Purity (GC) | ≥ 94.0 % |
| Inhibitor | 50 - 500 ppm |
| Phenothiazine | |

Specification: PRD.0.ZQ5.10000053309

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.