

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Methyl trans-vaccenate - ≥99% (capillary GC), liquid

Product Number:

V1381

CAS Number:

6198-58-9

MDL:

MFCD00064484

Formula:


C₁₉H₃₆O₂

Formula Weight:

296.49 g/mol

Storage Temperature:

-20 °C



TEST

Specification

Appearance (Turbidity)

Clear

Appearance (Color)

Colorless to Faint Yellow

Appearance (Form)

Liquid

Proton NMR spectrum

Conforms to Structure

Purity by Gas Chromatography

≥ 99 %

Specification: PRD.0.ZQ5.10000050557

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.