

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Methyl pyrimidine-5-carboxylate 95%

Product Number:

682241

CAS Number:

34253-01-5

MDL:


MFCD06204205

Formula:

C₆H₆N₂O₂

Formula Weight:

138.12 g/mol



TEST

Specification

Appearance (Color)	White to Orange to Brown
Appearance (Form)	Conforms to Requirements
Solid, Powder, Crystals, Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Purity (GC)	≥94.0 %
Carbon	49.0- 55.3 %
Nitrogen	19.1- 21.5 %

Remarks:

Specification Date : 11/03/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.