

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Biphenyl-3-carbonyl chloride - 97%

Product Number: **709913**
CAS Number: 42498-44-2

Formula: C₁₃H₉Cl
Formula Weight: 216.66 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Colorless to Yellow and White to Light Yellow
Appearance (Form)	Liquid or Solid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 96.5 %
Titration by AgNO ₃ After Oxygen Combustion	96.5 - 103.5 %

Specification: PRD.1.ZQ5.10000011089

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.