

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

3,4-Dimethoxystyrene - technical grade, contains 1% hydroquinone as inhibitor

Product Number:

154466

CAS Number:

6380-23-0

MDL:

MFCD00008615

Formula:


C₁₀H₁₂O₂

Formula Weight:

164.20 g/mol

Storage Temperature:

2 - 8 °C



TEST

Specification

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Liquid

Proton NMR spectra

Conforms to Structure

Purity (GC)

≥ 90.0 %

Inhibitor

0.25 - 1.75 %

1% hydroquinone as inhibitor

Specification: PRD.2.ZQ5.10000012894

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.