


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Biphenyl-4-sulfonyl chloride

Product Number: 536938
CAS Number: 1623-93-4
MDL: MFCD00078658
Formula: C₁₂H₉ClO₂S
Formula Weight: 252.72 g/mol



| TEST | Specification |
|----------------------------------|-----------------------|
| Appearance (Color) | White to Beige |
| Appearance (Form) | Powder |
| Infrared spectrum | Conforms to Structure |
| Purity (HPLC) | ≥ 90.0 % |
| Titration with AgNO ₃ | ≥ 90.0 % |

Specification: PRD.1.ZQ5.10000016009

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.