

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

Product Specification

Product Name: Monoclonal Anti-Cytokeratin, pan antibody produced in mouse – clone PCK-26, purified immunoglobulin, buffered aqueous solution  
Product Number: C5992

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	1.5 - 2.0 mg/ml
IHC - Tissue	Human placenta or skin
IHC - Fixation	FFPE-Formalin Fixed Paraffin Embedded
IHC - Working Concentration	10 - 20 µg/ml
IHC -Tissue	Human placenta or skin
IHC - Fixation	FFPE+ Antigen Retrieval by Decloaking
IHC - Working Concentration	10 - 20 µg/ml
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000017962

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.