


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Butyl methacrylate - 99%, contains monomethyl ether hydroquinone as inhibitor

Product Number: 235865
CAS Number: 97-88-1
MDL: MFCD00009444
Formula: C₈H₁₄O₂
Formula Weight: 142.20 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Refractive index at 20 °C	1.421 - 1.425
Purity (GC)	≥ 98.5 %
MEHQ as Inhibitor	8 - 28 ppm

Specification: PRD.2.ZQ5.10000023180

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.