


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Chloro-3-nitropyridine - 90%

Product Number: 674079
CAS Number: 13091-23-1
MDL: MFCD02683025
Formula: C₅H₃ClN₂O₂
Formula Weight: 158.54 g/mol
Storage Temperature: -20 °C



TEST

Specification

Appearance (Color)	White to Dark Yellow and White to Brown
Appearance (Form) Powder, Crystals or Crystalline Powder and/or Chunks	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Carbon	33.3 - 42.4 %
Nitrogen	15.5 - 19.8 %
Purity (HPLC)	≥ 88.0 %

Specification: PRD.1.ZQ5.10000010462

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.