

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Magnesium - in a Sure/Seal™ bottle, turnings, 37.5 mmol

Product Number: **723312**  
CAS Number: 7439-95-4

**Mg**

Formula: Mg  
Formula Weight: 24.31 g/mol

TEST	Specification
------	---------------

Appearance (Color)	Conforms to Requirements
Silver	
Appearance (Form)	Turnings
ICP Major Analysis	Confirmed
Confirms Magnesium Component	
Compleximetric EDTA	97.5 - 102.5 %

Specification: PRD.0.ZQ5.10000021269

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.