

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Phalloidin, Fluorescein Isothiocyanate Labeled - Phalloidin from Amanita phalloides

Product Number: P5282

MDL: MFCD00147902

Formula: C₅₈H₆₃N₁₀O₁₄S₄

Formula Weight: 1,252.44 g/mol

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Solubility (Color)

Faint Yellow to Dark Yellow and Faint
Green to Light Green and Yellow-Green
Clear

Solubility (Turbidity)

0.1 mg + 0.2 mL, MeOH

Composition by HPLC

mixed isomers

Conforms

Specification: PRD.1.ZQ5.10000036965

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.