


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Diethylene glycol - ReagentPlus®, 99%

Product Number: H26456
CAS Number: 111-46-6
MDL: MFCD00002882
Formula: C₄H₁₀O₃
Formula Weight: 106.12 g/mol



TEST	Specification
Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Liquid or Viscous Liquid
Infrared spectrum	Conforms to Structure
GC (area %)	≥ 98.5 %
Recommended Retest Period	-----
2 Years	
Note	Note
Product darkens over storage time without loss of purity	

Specification: PRD.3.ZQ5.10000031806

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.