

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Nickel on silica/alumina - extent of labeling: ~65 wt. % loading, powder

Product Number:

208779

Ni

MDL:

MFCD00011137

Formula:

Ni

TEST	Specification
------	---------------

Appearance (Color)

Black

Appearance (Form)

Powder

Gravimetric Analysis

61 - 69 %

% Ni, With Dimethylglyoxime

Specification: PRD.0.ZQ5.10000026247

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.