

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Halocarbon oil 27

Product Number: H8773

MDL: MFCD00147324



TEST	Specification
------	---------------

Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Viscosity	24 - 30
At 100 Degrees Fahrenheit	
Hydrocarbons	Pass
No More Than 10 PPM	
Saturation Test	Pass
No Less Than 15 Minutes	
Acidity	Pass
No More Than 0.0002 mg KOH/g	

Specification: PRD.0.ZQ5.10000021996

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.