


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
6-(p-Toluidino)-2-naphthalenesulfonic acid sodium salt

**Product Number:** T9792  
**CAS Number:** 53313-85-2  
**MDL:** MFCD00171528  
**Formula:** C<sub>17</sub>H<sub>14</sub>NNaO<sub>3</sub>S  
**Formula Weight:** 335.35 g/mol  
**Storage Temperature:** 2 - 8 °C



TEST	Specification
Appearance (Color)	White to Light Pink
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Solubility (Color)	Colorless to Light Yellow
Solubility (Turbidity)	Clear
10 mg + 4.0 ml; DMSO	
Purity (TLC)	≥ 98 %

Specification: PRD.0.ZQ5.10000039595

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.