


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Isosafrol - mixture of cis and trans

Product Number: 329606
CAS Number: 120-58-1
MDL: MFCD00005838
Formula: C₁₀H₁₀O₂
Formula Weight: 162.19 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
Refractive index at 20 ° C

Colorless to Yellow
Liquid
Conforms to Structure
1.574 - 1.578

Specification: PRD.1.ZQ5.1000005444

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.