SIGMA-ALDRICH®

sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name: Collagen from human placenta - Bornstein and Traub Type IV, powder

Product Number:	C8374
CAS Number:	9007-34-5

Storage Temperature: -20 ℃

TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear to Very Slightly Hazy
1 mg/mL, 0.5M HOAc	
Purity by SDS-PAGE	≥ 85 %
Greater than 85% of the collagen bands migrating with, or at a molecular weight greater than the alpha-1 band.	
Tested For Infectious Agents	Tested

Specification: PRD.0.ZQ5.10000041363

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.