


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,2-Dimethylsuccinic acid - 99%

Product Number: D186007
CAS Number: 597-43-3
MDL: MFCD00004196
Formula: C₆H₁₀O₄
Formula Weight: 146.14 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Crystal to Powder
Proton NMR Spectrum	Conforms to Structure
Carbon	48.6 - 50.0 %

Specification: PRD.0.ZQ5.10000044735

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.